

## Micro Positioning Systems with Piezo Inertial Drives

- Piezo inertial drive with lowest hysteresis
- Reached position will be held without current
- Long lifetime
- Variable step width
  - open loop about 200 nm
  - closed loop about 50 nm
- High pushing or blocking force
- No limit switches necessary
- Ultra compact dimensions possible
- Customized solutions as
  - ultra high vacuum uses
  - ultra low temperature uses
  - ultracompact drives
- Compact hand-held controller with
  - battery- or power supply operation
  - RS 232 Interface
- Compact USB controller with
  - USB voltage operation
  - USB - Interface



Miniature translation stage MS 30



Miniature translation stage MS 30; preparation for vacuum with linear measuring system



Customized development (positioning of an aperture in vacuum)



Monomode Coupler MK 25



Compact XYZ-positioner MX 25



Hand-held controller CN 30

# Table of contents of the delivery program

	<b>Page</b>
<b>1. Linear positioning systems with piezo inertial drives</b>	
1.1 Miniature Monomode Coupler MK 25	2
1.2 Linear positioners for one dimension	
1.2.1 Miniature Positioning Stage MS 15	4
1.2.1 Miniature Positioning Stage MS 30	6
1.2.2 Miniature Positioning Stage MS 38 with free aperture	8
1.2.3 Positioning Stage ML 17 for ultra low temperatures	10
1.3 Linear positioners for three dimensions	
1.3.1 Miniature Positioners MX 25 and MX 35	12
<b>2. Tilting systems with piezo inertial drives</b>	
2.1 Mirror tilting stage MT 25	14
<b>3. Position sensors for piezo inertial drives</b>	
3.1 Position Sensor PS 30 for Miniature Positioning Stage MS 30	16
<b>4. Motor Controllers</b>	
4.1 Controller for piezo inertial drives	
4.1.1 Handheld controller CN 30 for piezo inertial drives	18
4.1.2 USB – Controller CU 30 for piezo inertial drives	19
4.1.3 USB – Controller CF 30 for piezo inertial drives	20
4.2 Connecting cables for linear positioning systems	21
4.3 Autonomous controller CD 85 for DC-Motors	22
4.4 Autonomous controller CS 85 for stepper motors	25
<b>5. Customized Items</b>	
5.1 Linear measuring stage MS 30 with linear measuring system and vacuum	27
5.2 Motorized aperture for REM	27
5.3 XYZ – positioner for ultra low temperature microscope	27
5.4 Monomode positioner with MX 25	27

## MK 25

# Miniature Monomode Coupler MK 25 with piezo electric inertial motor



### Technical Data

Travel:	2 mm in xyz
Max. speed:	1.2 mm/s (depends on controller)
Optical height:	20 mm
Mass:	58 g

### Resolution (calculated)

Single step	~ 450 nm
1/16-step (with controller CU 30)	~ 30 nm
Half step	~ 250 nm
Double step (with controller CN 30)	~ 900 nm

### Lens type

Type A:	
Numerical aperture	0.40
Clear aperture	3.7 mm
Focal length	4.60 (670 nm) 4.64 (785 nm)
Type B:	
Numerical aperture	0.38
Clear aperture	2.5 mm

### Specifications

- Piezo driven step motor with low hysteresis
- holds reached position without current
- aspherical lenses with different NA
- optical height 20 mm
- 2 mm travel in xyz
- different fibre couplers
- step width about 200 nm
- positioning accuracy better than 100 nm with CU 30
- velocity up to 1.2 mm/s
- CNC-machined aluminium body
- no limit switches necessary
- vacuum preparation optionally
- customized designs possible
- driven by hand-held (CN.030.0001)  
or USB controller (CU.030.xx0x)

### Application Examples


- Micro-/Nano Technology
- Bio Technology
- Microscopy
- Telecommunication
- Metrology
- R & D

### Miniature Monomode- Coupler MK 25

Part no.


**MK.025.1**

	0	1	
Coupler type	FC	SMA	
	0	6	x
Vacuum	no	10 <sup>-6</sup>	customized
	A	B	x
Numerical aperture	0.40	0.38	customized

**Dimensions of the Monomode Coupler MK 25**

## MS 15

# Ultra small Miniature Translation Stages with piezo electric inertial drive



### Specifications

- Piezo driven step motor with low hysteresis
- holds reached position without current
- step width about 15 nm (depends on controller)
- positioning accuracy better than 50 nm
- velocity up to 1.5 mm/s (depends on controller)
- travel 3.5 mm
- xy or xyz combinations
- CNC-machined aluminium body
- precision ball bearing guides
- no limit switches necessary
- vacuum preparation optionally
- customized designs possible
- driven by hand-held (CN.030.0001) or USB controller (CU.030.xx0x) or USB controller (CF.030.xx0x)
- connection cable CK.030.1003 recommended

### Technical Data

Travel:	3.5 mm
Max. speed:	1.5 mm/s (depends on controller)
Mass:	6 g

### Load characteristics

Max. load	
$M_x, M_y, M_z$	0.4 Nm
$F_x$ (blocking force)	3 (4) N
$F_y, F_z$	20 N

### Resolution (calculated)

Single step	~ 450 nm
1/16-step (with controller CU 30)	~ 30 nm
1/64-step (with controller CF 30)	~ 10 nm
Half step	~ 250 nm
Double step (with controller CN 30)	~ 900 nm

### Guidance accuracy (without load)

Yaw angle	< 20 arc sec
Pitch angle	< 60 arc sec
Vertical deviation	< 1 $\mu$ m
Lateral deviation	< 2 $\mu$ m

### Application Examples

- Micro-/Nano Technology
- Bio Technology
- Microscopy
- Quality Control
- Metrology
- R & D

### Miniature Translation Stage MS 15

Part no.

**MS.015.04**

	0	6	9
Vacuum	no	$10^{-6}$	$10^{-9}$
	0	2	3
Mounting	x	xy (not removable)	xyz



**MS 30****Miniature Translation Stages  
with piezo electric inertial drive****Specifications**

- Piezo driven step motor with low hysteresis
- holds reached position without current
- step width about 10 nm (depends on controller)
- positioning accuracy better than 50 nm
- velocity up to 1.2 mm/s (depends on controller)
- travels up to 60 mm possible
- open- or closed loop-application
- xy or xyz combinations
- CNC-machined aluminium body
- precision ball bearing guides
- no limit switches necessary
- vacuum preparation optionally
- customized designs possible
- driven by hand-held (CN.030.0001)  
or USB controller (CU.030.xx0x)  
or USB controller (CF.030.xx0x)

**Application Examples**

- Micro-/Nano Technology
- Bio Technology
- Microscopy
- Quality Control
- Metrology
- R & D

**Technical Data**

Travel:	8, 18 or 30 mm
Max. speed:	1.2 mm/s (depends on controller)
Mass:	32, 38 or 54 g

**Load characteristics**

Max. load	
$M_x, M_y, M_z$	0.5 Nm
$F_x$ (blocking force)	4.5 (5) N
$F_y, F_z$	30 N

**Resolution (calculated)**

Single step	~ 450 nm
1/16-step	~ 30 nm
(with controller CU 30)	
1/64-step	~ 10 nm
(with controller CF 30)	
Half step	~ 250 nm
Double step	~ 900 nm
(with controller CN 30)	

**Guidance accuracy (without load)**

For 8 mm travel:	
Yaw angle	< 20 arc sec
Pitch angle	< 60 arc sec
Vertical deviation	< 1 $\mu$ m
Lateral deviation	< 2 $\mu$ m


**Miniature Translation  
Stage MS 30****Part no.****MS.030.**

	08	18	30
Travel	8 mm	18 mm	30 mm
	0	6	9
Vacuum	no	$10^{-6}$	$10^{-9}$
	0	1	
Measuring-system	without	with (please specify with sensor)	




## Series MS

## MS 30, 8 mm travel



## MS 30, 18 and 30 mm travel

Travel	A	B	C	D
18	40	33	10	15
30	52	45	16	21





## MS 38

# Miniature Translation Stages with free opening dia 10 mm and piezo electric inertial drive



### Specifications

- Piezo driven step motor with low hysteresis
- holds reached position without current
- step width about 15 nm (depends on controller)
- positioning accuracy better than 50 nm
- velocity up to 1.5 mm/s (depends on controller)
- travel up to 8 mm
- free opening of 10 mm (on whole travel)
- fits to optical bench systems (LINOS, OWIS aso.)
- xy or xyz combinations
- CNC-machined aluminium body
- precision ball bearing guides
- no limit switches necessary
- vacuum preparation optionally
- customized designs possible
- driven by hand-held (CN.030.0001) or USB controller (CU.030.xx0x)

### Technical Data

Travel:	8 mm
Free opening	10 mm
Max. speed:	1.5 mm/s (depends on controller)
Mass:	32 g

### Load characteristics

Max. load	
$M_x, M_y, M_z$	0.4 Nm
$F_x$ (blocking force)	4.5 (5) N
$F_y, F_z$	20 N

### Resolution (calculated)

Single step	~ 200 nm
1/16-step	~ 15 nm
(with controller CU 30)	

Half step	~ 100 nm
Double step	~ 400 nm
(with controller CN 30)	

### Guidance accuracy (without load)


Yaw angle	< 20 arc sec
Pitch angle	< 60 arc sec
Vertical deviation	< 1 $\mu$ m
Lateral deviation	< 2 $\mu$ m


### Application Examples

- Micro-/Nano Technology
- Bio Technology
- Microscopy
- Quality Control
- Metrology
- R & D

### Miniature Translation Stage MS 38

### Part no.

**MS.038.08□0**


**MS 38, 8 mm travel****xy and xyz combinations**

## ML 17

# Miniature Translation Stages for Ultra Low Temperatures with piezo electric inertial drive



### Technical Data

Travel:	5 mm
Max. speed:	1.0 mm/s (with controller CU 17 LT)
Electrical connection:	2 solder points
Mass:	25 g

### Load characteristics

Max. load	
$M_x$	3.0 Ncm
$M_y, M_z$	1.5 Ncm
$F_x$ (blocking force)	1.0 (1.5) N
$F_y, F_z$	1.0 N

### Resolution (calculated)

Single step	
with 22 V (at 4.2 K)	~ 200 nm
with 42 V (at 4.2 K)	~ 500 nm
with 82 V (at 4.2 K)	~ 1 $\mu$ m

(with controller CU 17 LT)

### Guidance accuracy (without load)


For 5 mm travel:	
Vertical deviation	< 2 $\mu$ m
Lateral deviation	< 2 $\mu$ m


### Specifications

- Piezo driven step motor with low hysteresis
- holds reached position without current
- for use at ultra low temperatures up to 4 K
- step width about 200 nm
- positioning accuracy better than 1  $\mu$ m
- velocity up to 1.0 mm/s
- travels up to 5 mm
- xy or xyz combinations possible  
(L-bracket ML.017.9001 for xyz needed)
- CNC-machined steel body
- precision linear bearings
- no limit switches necessary
- vacuum preparation optionally
- customized designs possible
- driven by USB controller (CU.017.0003)

### Application Examples

- Cryo - applications
- Micro-/Nano Technology
- Bio Technology
- Microscopy
- Quality Control
- Metrology
- R & D



**Dimensions of the ML 17****xyz – combination of ML 17**

## MX 25/35

# Miniature XYZ Positioners with piezo electric inertial drive



MX 25 with 2 mm travel in xyz-direction

## Specifications

- Piezo driven step motor with low hysteresis
- holds position without power
- up to 10 mm travel in xyz
- step width about 10 nm (depends on controller)
- positioning accuracy better than 50 nm
- velocity up to 1.2 mm/s (depends on controller)
- CNC-machined aluminium body
- no limit switches necessary
- vacuum preparation optionally
- customized designs possible
- driven by hand-held (CN.030.0001)  
or USB controller (CU.030.0003)  
or USB controller (CF.030.xx0x)

## Application Examples

- Micro-/Nano Technology
- Bio Technology
- Microscopy
- Quality Control
- Metrology
- R & D

## Technical Data

Travel:	
MX 25	2 mm in xyz
MX 35	10 mm in xyz
Max. speed:	1.2 mm/s (depends on the controller)
Mass:	
MX 25	46 g
MX 35	76 g

## Load characteristics

Type of load	
$M_x, M_y, M_z$	0.15 Nm
$F_x, F_y$ (jamming force)	< 2 N
$F_z$ (jamming force)	< 2 N

## Resolution (calculated)

Single step	~ 450 nm
1/16-step (with controller CU 30)	~ 30 nm
1/64-step (with controller CF 30)	~ 10 nm
Half step	~ 250 nm
Double step (with controller CN 30)	~ 900 nm



## Guidance accuracy (without load)

Vertical deviation	< 2 $\mu$ m
Lateral deviation	< 2 $\mu$ m

## Miniature XYZ Positioner Part no.

**MX.0 5.1 0**

	2	3	
Type	MX 25	MX 35	
	0	x	
Top plate design	Standard	custom design	
	0	6	9
Vacuum	no	10 <sup>-6</sup>	10 <sup>-9</sup>

**Drawings of the MX 25:****Drawings of the MX 35:**

## MT 25

# Kinematical Mirror Tilting Stage with piezo electric inertial drive



### Technical Data

Angular adjustment:	6 deg ( $\pm 3$ deg) in two orthogonal axes
Max. speed:	25 mrad/s (depends on controller)
Optical height:	11.5 mm (on rod) 12.5 mm (in plate)
Mirror mount:	dia. 12.7 mm (1/2 inch)
Free opening:	dia. 9 mm
Mass:	46 g
Max mass of mirror:	100 g

### Resolution (calculated)


Single step	$\sim 14 \mu\text{rad}$
1/16 step	$\sim 1 \mu\text{rad}$
(with controller CU 30)	
1/64-step	$\sim 0.3 \mu\text{rad}$
(with controller CF 30)	
Half step	$\sim 7 \mu\text{rad}$
Double step	$\sim 28 \mu\text{rad}$
(with controller CN 30)	

### Specifications




- Piezo driven step motor with low hysteresis
- holds reached position without current
- optical height 11.5 mm (on rod)
- angular adjustment 6 deg ( $\pm 3$  deg) in two axes
- for mirror with dia. 12.7 mm  
(with mirror adapter also for larger mirrors)
- free opening of 9 mm
- solid state hinges
- step width about  $0.3 \mu\text{rad}$  with CF 30
- velocity up to 25 mrad/s ( $> 1 \text{ deg/s}$ )
- also usable as prism stage
- no limit switches necessary
- vacuum preparation optionally
- customized designs possible
- driven by hand-held (CN.030.0001)  
or USB controller (CU.030.xx0x)  
or USB controller (CF.030.xx0x)

### Application Examples

- Micro-/Nano Technology
- Bio Technology
- Microscopy
- Telecommunication
- Metrology
- R & D





**Dimensions of the Mirror Tilting Stage MT 25****Coordinates of the two tilting axes****MT 25 with beamsplitter cube on optical rail**

## PS 30

# Linear Measuring System for Miniature Translation Stage MS 30



Miniature Translation Stage MS 30, 8 mm travel with PS 30

### Specifications

- measuring lengths up to 30 mm
- resolution better than 50 nm
- open- or closed loop-application
- xy combinations
- ultra high vacuum preparation optionally
- customized designs possible
- driven by USB-Controller CU 30 (CU.030.xx0x)

### Technical Data

Measuring length	8, 18 or 30 mm
Mass	about 20 g

### Resolution

Standard	1 $\mu\text{m}$
Premium	0.5 $\mu\text{m}$
High End	0.1 $\mu\text{m}$
Excellence	0.05 $\mu\text{m}$
(other resolutions on request)	

### Scale tape

Material	Steel
Grating period	20 $\mu\text{m}$
Length of scale	measuring length + 22 mm
Reference mark	at the center of the scale
Linear expansion coefficient	$23.8 \times 10^{-6} \text{ grd}^{-1}$
Accuracy class	$\pm 1 \mu\text{m/m}$

### Electrical Data

Scanning frequency	max. 400 kHz
Output signal	RS 422 with interpolation
Supply voltage	5 V DC $\pm 10\%$
Power consumption	150 mA
Cable length	1.5 m
Operating temperature	$0^\circ - 55^\circ \text{C}$
Vibration (50 – 2,000 Hz)	$< 200 \text{ ms}^{-2}$
Shock (11 ms)	$< 400 \text{ ms}^{-2}$


### Application Examples

- Micro-/Nano Technology
- Bio Technology
- Microscopy
- Quality Control
- Metrology
- R & D

### Linear Measuring System Part no. PS30 for Miniature Stage MS 30

PS.030.□□□□

	08	18	30	
Travel	8 mm	18 mm	30 mm	
	0	3	6	9
Vacuum	no	$10^{-3}$	$10^{-6}$	$10^{-9}$
	2	3	5	6
Resolution	1 $\mu\text{m}$	0.5 $\mu\text{m}$	0.1 $\mu\text{m}$	0.05 $\mu\text{m}$

**Linear Measuring System PS 30 with adapter plate for MS 30****XYZ-arrangement with MS 30 (8 mm travel) and PS 30 (8 mm travel and 50 nm resolution)**

**CN 30****Controller for Piezo Electric Inertial Motors****Specifications**

- Special design for piezo driven step motor
- works with batteries or external power supply
- operates up to 3 axes
- customized designs possible
- USB-interface optionally

**Technical Data**

- Power supply:** 4 x round cell LR6 (Alkaline ~ 2000 mAh) and portable power supply  
in: 100 – 240 V  
out: 6 VDC with user-plug  
Ø 5.5/2.1 mm
- Operation modes:** "Local" with keyboard (current consumption ~ 20 µA)  
"Remote" with RS 232-Interface (current consumption ~6 mA incl. LED)  
change switch for both modes on keyboard
- Connections:** 9pin Sub D-connector (male)  
9pin Sub D-connector (female)  
2pin Socket (female)
- Speed modes:** "Fast" (ca. 0.7 mm/s)  
"Slow" ("Fast"/2)  
"Single step"
- Step modes:** "Full step"  
"Half step" ("full step"/2)  
"Double step" ("full step" x 2)
- RS 232-interface to the internal controller**
- Connection:** 9pin Sub D-connector (male)  
**Data rate:** 19200 Baud
- Software for PC**
- for DOS, Win 95, Win 2000, Win XP and Win NT: executable software with action buttons (similar to the stand-alone keyboard)
  - for use with own programs: Borland Pascal Unit (DOS)  
DLL for Windows
- Mass:** 370 g incl. batteries
- Dimensions:** 158 x 96 x 52 mm
- Delivery includes:** - miniature power supply  
- 4 x round cell LR6  
- RS 232 connecting cable

**Application Examples**

- Micro-/Nano Technology
- Bio Technology
- Microscopy
- Quality Control
- Metrology
- R & D

**CN 30 Controller****Part no.****CN.030.0001**

## CU 30

# USB-Controller for Piezo Electric Inertial Drives



CU 30 without sensor (for three axes)

### Specifications

- Special design for piezo driven inertial motors
- works with USB-voltage (without sensor) or an additional power supply (closed loop controller)
- applications without position sensor or closed loop
- operates up to 3 axes
- customized designs possible (interfaces or number of axes)




CU 30 closed loop (for three axes)

### Application Examples

- Micro-/Nano Technology
- Bio Technology
- Microscopy
- Quality Control
- Metrology
- R & D

### Technical Data

Power supply:	USB-interface (without sensor) or portable power supply (for CU 30CL) in: 100 – 240 V out: 6 VDC with user-plug
Operation:	via USB-Interface 1.1 (without sensor) via USB-Interface 2.0 (closed loop)
Connections:	
without sensor	1 x 9pin Sub D-connector (female) 1 x USB Type B socket (USB 1.1)
closed loop	3 x 2pin MMCX sockets (female) 3 x 9pin Sub D-connectors (female) 1 x USB Type B socket (USB 2.0)
Speed modes:	0 to 1000 (0 to ca. 1.2 mm/s) for CU 30 0 to 255 (0 to ca. 1.2 mm/s) for CU 30CL
Software for PC	
- for DOS, Win 95, Win 2000, Win XP and Win NT:	executable software with action buttons
- for use with own programs:	Borland Delphi 7.0 Unit DLL for Windows Wrapper-DLL for LabView etc.
Mass:	about 150 g (without sensor) (also closed loop for 1 axis) about 470 g (closed loop for 3 axes)
Dimensions:	
without sensor	118 x 86 x 26 mm (L x W x H)
closed loop	124 x 170 x 55 mm (L x W x H)
Delivery includes:	USB connecting cable (CK.030.USB0) power supply (for 3 axes closed loop)



## CF 30

# USB-Controller for Piezo Electric Inertial Drives



CF 30 (for three axes)

### Specifications


- Special design for piezo driven inertial motors
- Works with USB-voltage only (one axis move at the same time)
- Works with an additional power supply (simultaneous move off three axes possible)
- Auto detection of the power supply
- Operates up to 3 axes
- Plain text ASCII commands for easy integration in your own environment
- Drivers for many operating systems
- Every full step has 64 micro steps
- Special wave forms for smooth movement
- customized designs possible (interfaces or number of axes)

### Technical Data

Power supply:	USB-interface (one axis) or portable power supply in: 100 – 240 V out: 6 VDC with user-plug
Interface (host):	Via USB 1.1 Via USB 2.0 (Full speed)
Interface (client):	USB to serial converter FT232R from FTDI ( <a href="http://www.ftdichip.com">www.ftdichip.com</a> )
Connections:	Stage: 9pin Sub D-connector (female) PC: USB Type B socket (USB 1.1)
Speed modes:	-32767 to 32767 (0 to about 1.5 mm/s) Sawtooth voltage (0 to 3.1 KHz) Up to 800 commands/s
Software for PC	
- for	Win XP and Vista (32 bit): Executable software with GUI
- for use with own programs:	Drivers from FTDI: Future Technology Devices International Ltd. Plain text ASCII commands FTDI Drivers for Windows, Linux, Mac OS X, ....
Mass:	about 195 g
Dimensions:	118 x 86 x 26 mm (L x W x H)
Delivery includes:	USB connecting cable (CK.030.USB0) Power supply (for 3 axes)

### Application Examples

- Micro-/Nano Technology
- Bio Technology
- Microscopy
- Quality Control
- Metrology
- R & D



**CK 25/30****Connecting Cables for Piezo Electric Inertial Drives**


Part-no.	Description	Connections	Length
CK.025.0002	connecting cable for 1 x MT 25 (standard)	9pin Sub D-connector (male) to 3pin inline plug	1 m
CK.025.0003	connecting cable for 1 x MK 25 or 1 x MX 25/35 (standard)	9pin Sub D-connector (male) to 4pin inline plug	1 m
CK.025.1002	connecting cable for 1 x MT 25 (open wires)	9pin Sub D-connector (male) to 3 open wires	1 m
CK.025.1003	connecting cable for 1 x MK 25 or 1 x MX 25/35 (open wires)	9pin Sub D-connector (male) to 4 open wires	1 m
CK.025.1062	vacuum connecting cable (up to $10^{-6}$ mbar) for 1 x MT 25 (open wires)	3pin inline plug to to 3pin inline plug	1 m
CK.025.1063	vacuum connecting cable (up to $10^{-6}$ mbar) for 1 x MK 25 or 1 x MX 25/35 (open wires)	4pin inline plug to 4pin inline plug	1 m
CK.025.1092	vacuum connecting cable (up to $10^{-9}$ mbar) for 1 x MT 25 (open wires)	3pin inline plug to directly mounting to MT 25	1 m
CK.025.1093	vacuum connecting cable (up to $10^{-9}$ mbar) for 1 x MK 25 or 1 x MX 25/35 (open wires)	4pin inline plug to directly mounting to MX 25/35	1 m
CK.030.0001	connecting cable for 1 x MS 30 (for closed loop controller CU 30 CL)	1 MMCX-plug (angled) to 1 MMCX-plug (angled)	1 m
CK.030.0001.3m	connecting cable for 1 x MS 30 (for closed loop controller CU 30 CL)	1 MMCX-plug (angled) to 1 MMCX-plug (angled)	3 m
CK.030.0003	connecting cable for up to 3 x MS 30 (standard)	9pin Sub D-connector (male) to 3 MMCX-plugs (angled)	1 m
CK.030.0003.2m	connecting cable for up to 3 x MS 30 (standard)	9pin Sub D-connector (male) to 3 MMCX-plugs (angled)	2 m
CK.030.0003.3m	connecting cable for up to 3 x MS 30 (standard)	9pin Sub D-connector (male) to 3 MMCX-plugs (angled)	3 m
CK.030.1061	vacuum connecting cable (up to $10^{-6}$ mbar) for 1 x MS 30 (open wires)	2pin inline plug to to 1 MMCX-plug	1 m
CK.030.1091	vacuum connecting cable (up to $10^{-9}$ mbar) for 1 x MS 30 (open wires)	2pin inline plug to directly mounting to MS 30	1 m
CK.030.USB0	RS 232 – USB Transformer (only applicable with CN 30)	9pin Sub D-connector (male) to USB A – plug	1 m
CK.030.USB0	USB connecting cable	USB A – to USB B connectors	1.8 m

For vacuum connections please ask with detailed technical specifications.



## CD 85

## Controller for Servo Motors



### Specifications


- compact design, only 146 x 85 x 55 mm
- for up to 4 motors with 100 Watts
- power supply 24 Volts
- I/O card with 16 channels optionally


### Technical Data

No. of motors:	1 to 4 motors
Types of motors:	brushed or brushless servo
Encoder:	incremental
Types of profiles:	trapezoidal, s-curve, velocity profile, electrical gear, closed curves
Contr. commands:	user friendly and self-instructional commands and parameters
Electronical gear:	master encoder selectable, slave-axis selectable, factor "in motion" alterable
Outline buffer:	any outline with up to 65,000 points programmable
Control parameter:	PID-portion, acceleration/velocity parameter changeable during motion
Monitoring functions:	error limit and position error window selectable, comparator function, simultaneous tracing of different parameters during motion
Inputs per axis:	Encoder signals A, B, Index, Home, End+, End-
Free inputs:	16 inputs, protected against misconnection, short-circuit proof
Free outputs:	16 outputs
Status/Error display:	6 LED
Interfaces:	RS 232, CAN-Bus (standard CAN 2.0B)
Supply voltage:	24 V

### Application Examples


- Motion control in
- Bio Technology
  - Microscopy
  - Quality Control
  - Metrology
  - R & D



**Dimensions:**

## CS 85

## Controller for Stepper Motors



### Technical Data

No. of motors:	1 to 4 motors
Types of motors:	stepper or micro stepper motors
Encoder:	incremental
Types of profiles:	trapezoidal, s-curve, velocity profile, electrical gear, closed curves
Contr. commands:	user friendly and self-instructional commands and parameters
Electronical gear:	master encoder selectable, slave-axis selectable, factor "in motion" alterable
Outline buffer:	any outline with up to 65,000 points programmable
Control parameter:	PID-portion, acceleration/velocity parameter changeable during motion
Monitoring functions:	error limit and position error window selectable, comparator function, simultaneous tracing of different parameters during motion
Inputs per axis:	Encoder signals A, B, Index, Home, End+, End-
Free inputs:	16 inputs, protected against misconnection, short-circuit proof
Free outputs:	16 outputs
Status/Error display:	6 LED
Interfaces:	RS 232, CAN-Bus (standard CAN 2.0B)
Supply voltage:	24 V

### Specifications


- compact design, only 146 x 85 x 55 mm
- for up to 4 motors with 100 Watts
- power supply 24 Volts
- I/O card with 16 channels optionally


### Application Examples

- Motion control in
- Bio Technology
  - Microscopy
  - Quality Control
  - Metrology
  - R & D

### Controller for Stepper Motors Part no.


CS.085.□00□



**Dimensions:**

## Customized items of the mechOnics ag

### Linear measuring stage MS 30 with linear measuring system and vacuum preparation




MS 30 with

- 8 mm travel
- vacuum preparation for  $10^{-6}$  mbar

Linear measuring system PS 30 with

- grating period  $20\text{ }\mu\text{m}$
- measuring increments  $0.1\text{ }\mu\text{m}$
- vacuum preparation for  $10^{-6}$  mbar

### Motorized aperture for REM



Linear positioner for an aperture with


- 7 mm travel in x – direction
- 1 mm travel in y – direction
- positioning accuracy better  $1\text{ }\mu\text{m}$
- vacuum preparation for  $10^{-6}$  mbar
- nonmagnetic materials

### XYZ – Positioner for ultra low temperature microscope



- 3 mm travel in xyz
- About 250 nm step width
- Working temperature 4.2K
- vacuum preparation for  $10^{-6}$  mbar
- Compact design

### XYZ – Positioner MX 25 as Fiberpositioner



XYZ – Positioner MX 25 with

- 2 mm travel
- About 200 nm step width
- Velocity about 0.5 mm
- On profile rail with optical height 40 mm